

PEN-PRO DIGITAL REFRACTOMETER

EASY BRIX MEASUREMENT WITH THIS HAND HELD 'PEN' REFRACTOMETER BY 



- Water Resistant (IP65)
- Auto Temperature Compensation
- Calibration With Water
- Measurement in 2 Seconds
- Light and Compact (70g)

The PEN-PRO has a measurement range of Brix 0.0 to 85.0% which can be used to measure almost any sample.

There are 2 ways to measure with the PEN-PRO. Either dip the tip into the sample and press the START key OR press the START key and touch the tip into the sample.

There is also the new feature to perform continuous measurements.

The temperature range of the "Dip-Type" Refractometer PEN-PRO has been improved from 60°C to 100°C. No more wait ! Measure hot samples while in the cooking process, such as soups and sauces.



Easy continuous measurements

Once the START key is pressed, measurements are taken every second.

Improved efficiency



Measure different samples one after another, saving time and improving efficiency.

Easy to clean



The PEN-Pro is fully waterproof and can be washed under the tap

- Scale Range: 0%-85% Brix
- Resolution: 0.1% Brix
- Accuracy: +/- 0.2%
- Dimensions: 16cm x 3.8cm x 1.8cm
- Weight: 70 grams/3 ounces
- Sample Temperature: 10 to 100°C
- Power Supply: 1 x AAA Battery

LABCELL LTD

FOUR MARKS, ALTON, HAMPSHIRE GU34 5PZ
TEL: ++44 (0)1420 568150 FAX: ++44 (0)1420 568151
e: mail@labcell.com www.labcell.com



Scan the QR Code with your Smartphone to learn more about the Digital Refractometers